



V&Q MANUFACTURING CORPORATION

Add: 4A, Fuyuan Road, Pengchang town, Xiantao City, Hubei Province, China
Tel: 0086 27 87680925 Fax: 0086 27 87680926 Web: www.chinavqmc.com

Data Sheet For Face Mask



1. Item : Medical FACE MASK (Type II R)

2. Picture:



3. Requirements for EN14683:2019:

Table 1 — Performance requirements for medical face masks

Test	Type I ^a	Type II	Type IIR
Bacterial filtration efficiency (BFE), (%)	≥ 95	≥ 98	≥ 98
Differential pressure (Pa/cm ²)	< 40	< 40	< 60
Splash resistance pressure (kPa)	Not required	Not required	≥ 16,0
Microbial cleanliness (cfu/g)	≤ 30	≤ 30	≤ 30

^a Type I medical face masks should only be used for patients and other persons to reduce the risk of spread of infections particularly in epidemic or pandemic situations. Type I masks are not intended for use by healthcare professionals in an operating room or in other medical settings with similar requirements.



4. Test Reports:

4.1 Bacterial Filtration Efficiency (B.F.E.) : 99.4% average

(Tested by Nelson Laboratories Inc.)

4.2 Particle Filtration Efficiency (P.F.E.) : 99.74% average

(Tested by Nelson Laboratories Inc.)

4.3 Differential Pressure (Delta P) : 2.7 mm-H₂O average.

(Tested by Nelson Laboratories Inc.)

4.4 Synthetic Blood Penetration: Pass at 120 mm Hg

(Tested by Nelson Laboratories Inc.)

4.5 Flammability Test : Pass

(Tested by Nelson Laboratories Inc.)

4.6 Skin irritation Test: No irritation to skin

(Tested by Soochow University)

4.7 Skin Sensitization Test: No Sensitization to skin

(Tested by Soochow University)

4.8 In Vitro Cytotoxicity Test: No Cytotoxicity

(Tested by Soochow University)

5. Size : 3-ply, 17.5x9.5 cm

6. Material :

Outer layer: polypropylene

Middle layer: Melt blown

Inner layer: polypropylene

7. Nose Bridge :

11cm (long) X 0.2cm (wide) Adjustable Aluminum-plastic Band.

8. Elastic Band : latex-free elastic band.

9. Color :Blue

10. Packing Details :

50 pieces of masks in a box.

20 boxes in a carton.





V&Q MANUFACTURING CORPORATION

Add: 4A, Fuyuan Road, Pengchang town, Xiantao City, Hubei Province, China
Tel: 0086 27 87680925 Fax: 0086 27 87680926 Web: www.chinavqmc.com

11. Shelf Life: 5 years from the date of production

12. Instructions for use

- 1). Check the face mask before using. It should be non-polluting, non-damaging, ear loop with good elasticity and nose bridge without break.
- 2). It is disposable product, and shouldn't be used recycling.
- 3). It should not be used more than 12 hours.
- 4). Don't touch the inner side during using.
- 5). Don't use if any breakage of packing or product expired.
- 6). Please dispose the mask as medical waste once used.

13. Handling and Storage

Don't store in conditions of dampness or high temperature.

MINIMAL TEMPERATURE °C is 10 °C

MAXIMUM TEMPERATURE °C is 35 °C


MINIMAL RELATIVE HUMIDITY (%) is 45%

MAXIMUM RELATIVE HUMIDITY (%) is 65%

14. EU Representative

Caretechion GmbH

Address: Niederrheinstr. 71, 40474 Duesseldorf, Deutschland



Jacky Zou