



2024

Catalogue

Single-use medical devices

APPLIED SA

INDEX

Sterile Product

Non sterile Product


PRESENTATION	
List of symbols	3
About the company	4
History	4
"Tailor-made" set	5
The advantages of single use	6

MEDICAL DEVICES & ACCESSORIES	
Personal protective equipment	
<i>Gowns - overalls</i>	60
<i>Anti-slip shoe cover, cap & protective clothing set</i>	61
Protective bed clothes	62

MEDICAL DEVICES & ACCESSORIES	
Applicators	
<i>Applicators wood</i>	7-8
<i>Applicators polypropylene (PP)</i>	9-10
<i>Applicators polystyrene (PS)</i>	11-14
<i>Applicators aluminium</i>	15-16
<i>Applicators stainless steel</i>	17-18
Copolymer gloves	19
Compresses	
<i>Woven compresses</i>	20-23
<i>Non-woven compresses</i>	24-25
<i>Pre-cut non-woven compresses</i>	26
<i>Abdominal compresses</i>	27
Gauze balls	
<i>Gauze balls</i>	28
<i>Gauze & non-woven balls</i>	29
<i>Miscellaneous balls</i>	30
<i>Mini tampons</i>	31
Tongue depressor	32
Kidney-shaped bowl	33-34
Surgical drapes	35
Cups & containers	36
Plastic instruments	37
Scalpels	38
Metal instruments	
<i>Scissors</i>	39-40
<i>Hook & curettes</i>	41
<i>Retractors & speculum</i>	42
<i>Forceps</i>	43-44
<i>Tweezers</i>	45-46
<i>Needle holders & blades</i>	47
<i>Gynaecology</i>	48-49
<i>Washing process of metal instruments</i>	50
<i>Bipolar forceps</i>	51-53
<i>Cannulas</i>	
<i>Olive tip - stainless steel</i>	54-55
<i>Blunt tip - stainless steel</i>	56
<i>Olive tip - plastic</i>	57
Collection and recycling of single-use metal instruments	58
Surgical face masks	59

SETS	
Dressing sets	63-68
Disinfection sets	69-73
Skin disinfection/preparation set	74
Wound care dressing sets	75-87
Intravenous start set	88-90
Suture sets	91-102
Suture removal sets	103-105
Catheterisation sets	106-109
Minor surgery dermatology set	110
Basic surgery set	111
Surgical set	112
Minor facial surgery set	113
Hand surgery tray sets	114-115
Phlebectomy set	116
Ingrown nail sets	117-118
Maternity suture set	119
Episiotomy set	120
IUD insertion set	121
Endoscopy set	122-123
Endoscope cleaning care sets	124-125
CT scanner injection set	126
Precision intervention set	127
Intervention set	128
Abscess set (standard)	129
Basic set	130
Biopsy Set	131
Diabetology set	132
Drainage set	133
Injection set (standard)	134
Multipurpose set	135
Puncture set	136
Small puncture set	137
Emergency set	138

The list of products referred above is not exhaustive.



LIST OF SYMBOLS



Sterilised using ethylene oxide



Non-sterile



Do not reuse
Single use



Do not use if packaging is damaged



European conformity marking



CE certification issued by the notified body with accreditation ID 1639 (SGS - Belgium) for marking under Directive 93/42/EEC



Distributed product



Refer to the instructions for use



1st level of packaging (pieces, simple packaging, peelpack)



2nd level of packaging (box)



3rd level of packaging (carton)

ABOUT THE COMPANY

The Applied SA Group, **manufacturer of sterile and non-sterile single-use medical devices**, is known primarily as one of the largest manufacturers of medical applicators in Europe.

We are certified EN ISO 13485:2016 and our medical devices carry the CE 1639 marking under Annex V of the Medical Devices Directive 93/42/EEC for sterility (class I-s) and under Annex II, section 4 for medical devices in class IIa. Applied SA's main customers include distributors of medical devices, hospitals and clinics, care/nursing homes for elderly, private medical practices, home care services, and pharmacies.




Our know-how, experience and flexibility are all assets that make our company a leader in its field. We can offer you a wide variety of single-use medical devices and various associated services such as ethylene oxide sterilisation assured by the C.S.S.R (Centre de Stérilisation de Suisse Romande), certified EN ISO 13485.

Applied SA makes a major contribution to the comfort and safety of patients and healthcare professionals through its work on producing and bringing to the market single-use medical devices which are perfectly suited to the needs of the medical sector.




HISTORY





«TAILOR-MADE» SET

Applied SA, with its extensive medical knowledge, manages the entire process of creating your single use custom set from start to finish, ensuring that the final product satisfies all of your needs in terms of quality of care, safety and management constraints.




Our specialists are available to provide any information you require or to discuss any project you may wish to develop.

The development and production of tailor-made sets are subject to minimum quantities.




THE ADVANTAGES OF SINGLE USE



- Single use is a safer option because cross-contamination is avoided.
- Single use eliminates the problem of cleaning reusable equipment and reduces costs associated with logistics (it saves space), execution (it saves considerable time), maintenance and investment.
- Choosing single-use products reduces the carbon footprint by saving thousands of litres of water and detergents. Since we are concerned about the impact on the environment, we provide our customers with a free service for the recovery and recycling of our metal instruments after use.

**EXPERTISE AND EXPERIENCE OF PRODUCING
«SINGLE-USE MEDICAL DEVICES»**




CE
1639

STERILE EO




1102245

Applicator wood
150x2.2mm
Small head 5-6mm
Cotton




1102246

Applicator wood
150x2.2mm
Small head 5-6mm
Cotton




1101274

Applicator wood
150x2.2mm
Small head 5-6mm
Cotton




1101277

Applicator wood
200x2.5mm
Small head 5-6mm
Cotton




1111413

Applicator wood
150x2.5mm
Medium head 7-8mm
Cotton




1121243

Applicator wood
150x3mm
Big head 11-12mm
Continental shape, cotton




1121463

Applicator wood
150x3mm
Big head 11-12mm
American shape, cotton




1121469

Applicator wood
150x3mm
Big head 11-12mm
Continental shape, cotton




1121475

Applicator wood
150x3mm
Big head 11-12mm
American shape, cotton




1131563



Applicator wood
200x3mm
Big head 11-12mm
Continental shape, cotton



1131461


Applicator wood
200x3mm
Mega head 13-15mm
American shape, cotton






1001202

Applicator wood
150x2.2mm
Small head 5-6mm
Cotton




1001230

Applicator wood
150x2.2mm
Small head 5-6mm
Cotton




1011124

Applicator wood
150x2.5mm
Medium head 7-8mm
Cotton




1021275

Applicator wood
150x2.5mm
Big head 11-12mm
Continental shape, cotton




1021292


Applicator wood
150x3mm
Big head 11-12mm
American shape, cotton



1022467

Applicator wood
200x3mm
Big head 11-12mm
Continental shape, cotton





CE
1639


STERILE EO



1323000

Applicator PP
75x2.5mm
Small head 5-6mm


Viscose



1501256


Applicator PP
150x2.5mm
Small head 5-6mm

Cotton




1521288

Applicator PP
150x4mm
Big head 11-12mm
American shape, cotton




1521289

Applicator PP
150x4mm
Big head 11-12mm
Continental shape, cotton




1521296

Applicator PP
150x4mm
Big head 11-12mm
Continental shape, cotton




1522252



Applicator PP
150x4mm
Big head 11-12mm
Continental shape, cotton



1532290


Applicator PP
200x4mm
Big head 11-12mm
Continental shape, cotton






1401219

Applicator PP
150x2.5mm
Small head 5-6mm
Cotton




1402244

Applicator PP
150x2.5mm
Small head 5-6mm
Viscose




1421283

Applicator PP
150x4mm
Big head 11-12mm
American shape, cotton




1421504

Applicator PP
300x4mm
Big head 11-12mm
Continental shape, cotton




1422278

Applicator PP
150x4mm
Big head 11-12mm
Continental shape, cotton




1422280

Applicator PP
150x4mm
Big head 11-12mm
American shape, cotton




1422286

Applicator PP
100x4mm
Big head 11-12mm
Continental shape, cotton




1422287

Applicator PP
100x4mm
Big head 11-12mm
American shape, cotton




1422505

Applicator PP
200x4mm
Big head 11-12mm
Continental shape, cotton



1432504

Applicator PP
397x4mm
Mega head 13-15mm
American shape, cotton





APPLIMED SA

APPLICATORS Polystyrene (PS)




CE
1639

STERILE EO




1301241

Applicator PS
150x2.5mm
Small head 5-6mm
Viscose



1301265


Applicator PS
150x2.5mm
Small head 5-6mm
Polyester






APPLIED SA

APPLICATORS Polystyrene (PS)




1211246

Applicator PS
150x2.5mm
Small head 5-6mm
Polyester, semi choc




1202180

Applicator PS
150x2.5mm
Small head 5-6mm
Cotton, semi choc




1202218

Applicator PS
150x2.2mm
Small head 5-6mm
Viscose




1211247

Applicator PS
93x2.5mm
Small head 5-6mm
Viscose



1211248


Applicator PS
125x2.5mm
Small head 5-6mm
Viscose, semi choc





APPLIMED SA

APPLICATORS
Polystyrene (PS)




CE
1639

STERILE EO



1251152


Applicator PS
140x4mm
Green foam - **Choc**






APPLIED SA



APPLICATORS Polystyrene (PS)



1251153

Applicator PS
140x4mm
Green foam - **Choc**






CE
1639

STERILE EO




1901199


Applicator aluminium
150x0.9mm
Small head 2-3mm
Cotton



1901941


Applicator aluminium
150x0.9mm
Small head 2-3mm
Polyester






1801183

Applicator aluminium
150x0.9mm
Small head 2-3mm
Cotton




1801193



Applicator aluminium
150x0.9mm
Small head 2-3mm
Viscose



1806197

Applicator aluminium
150x0.9mm
Small head 2-3mm
Polyester





CE
1639

STERILE EO




2101188

Applicator stainless steel

150x0.9mm

Mini head 2.2-2.7mm


Cotton





APPLIMED SA

APPLICATORS
Stainless steel



CE




2001192


Applicator stainless steel

150x0.9mm

Mini head 2.2-2.7mm

Cotton






STERILE	EO
---------	----




3109114


Copolymer gloves
Size M



3109115

Copolymer gloves
Size L






CE
1639

STERILE EO




4501992

12-ply compress
10x10cm, woven,
without X-ray thread




4502037

8-ply compress
10x10cm, woven,
without X-ray thread




4502043

16-ply compress,
10x10cm, woven,
without X-ray thread
Double packaging




4502047

12-ply compress,
5x5cm, woven,
without X-ray thread
Double packaging




4502065

12-ply compress,
5x5cm, woven,
without X-ray thread




4502075

16-ply compress,
10x10cm, woven,
with X-ray thread
Double packaging




4502079

12-ply compress,
10x10cm, woven,
without X-ray thread
Double packaging




4801126


8-ply compress,
1x normal
2x 25mm incision 5x5cm
woven,
without X-ray thread



4901893


8-ply compress,
10x10cm, woven,
without X-ray thread






4502029

8-ply compress,
10x10cm, woven,
without X-ray thread




4502051

12-ply compress,
5x5cm
without X-ray thread




4502052


12-ply compress,
10x10cm
without X-ray thread



4502083

8-ply, compress,
5x5cm
without X-ray thread






CE
1639

STERILE EO




4501990

24-ply compress,
10x20cm, woven,
without X-ray thread




4502038

24-ply compress,
10x20cm, woven,
without X-ray thread
Double packaging




4502039

24-ply compress,
10x20cm, woven,
without X-ray thread
Double packaging




4502041


24-ply compress,
10x20cm, woven,
with X-ray thread
Double packaging



4801127

12-ply compress,
7.5x10cm, woven,
«I» incision






APPLIMED SA


Woven COMPRESSES



4502053

12-ply compress,
10x20cm, woven,
without X-ray thread






STERILE EO




4801145

4-ply compress
5x5cm




4801146

4-ply compress
7.5x7.5cm




4801147

4-ply compress
10x10cm




4801148

4-ply compress
10x20cm




4801809

4-ply compress
5x5cm



4802206


4-ply compress
10x10cm
Double packaging




CE
1639

4802208


4-ply compress
5x5cm
Double packaging






4702111

4-ply compress
5x5cm




4702112

4-ply compress
7.5x7.5cm




4702113


4-ply compress
10x10cm



4702115

4-ply compress
10x20cm





APPLIMED SA

Pre-cut non-woven COMPRESSES




CE
1639


STERILE EO



4801132

4-ply, compress
5x5cm
Pre-cut «Y»





APPLIMED SA

Abdominal COMPRESSES



CE
1639

STERILE EO




4502059


Compress

6 layers, 45x45cm

**Prewashed with lace
and X-ray thread**


Double packaging






4002030

Gauze balls
Ø 3cm
without X-ray thread




4002035

Gauze balls
Ø 3.5cm
without X-ray thread



4002050


Gauze balls
Ø 5cm
without X-ray thread



4002425

Gauze balls
Ø 2.5cm
with X-ray thread






CE
1639

STERILE EO




4102002

Gauze balls
Ø 2.5 - 3 cm
without X-ray thread




4102008

Gauze balls
Ø 2.5 - 3 cm
without X-ray thread




4102009

Gauze balls
Ø 7 - 8 cm
without X-ray thread




4102263

Gauze balls
Ø 4 - 4.5 cm
with X-ray thread
Double packaging




4102264

Gauze balls
Ø 4 - 4.5 cm
without X-ray thread
Double packaging




4102265

Non-woven balls
Ø 4.5 cm
without X-ray thread
Double packaging




4102266

Gauze balls
Ø 3 cm
with X-ray thread
Double packaging




4102313

Non-woven balls
Ø 3 - 4 cm
without X-ray thread
Double packaging




4802118


Gauze balls
Ø 3 cm
with X-ray thread



4802212

Gauze balls
Ø 4 - 4.5 cm
without X-ray thread





CE
1639

STERILE	EO
---------	----




4101100


Non-woven balls

Ø 4.5 cm

without X-ray thread

In a rigid blister pack with
2 compartments






CE
1639

STERILE EO




4101105

Woven mini tampon
Micro-surgery
10x10cm
with X-ray thread RX
Double packaging




4101140

Woven mini tampon
Micro-surgery
4x4cm
with X-ray thread RX
Double packaging




4102160

Woven mini tampon
Micro-surgery
10x10cm
with X-ray thread RX
Double packaging




4102161

Woven mini tampon
Micro-surgery
8x8cm
with X-ray thread RX
Double packaging




4102310

Woven mini tampon
Micro-surgery
4x4cm
with X-ray thread RX
Double packaging



4102311

Woven mini tampon
Micro-surgery
6x6cm
with X-ray thread RX
Double packaging





APPLIMED SA

TONGUE DEPRESSOR




CE
1639

STERILE	EO
---------	----



5101802

Wooden tongue depressor
1.8x15cm






STERILE	EO
---------	----



4801019


Kidney-shaped bowl, sterile
APET, 252 x 150 x 40 mm







9504000

Kidney-shaped bowl,
APET, 252 x 150 x 40 mm






CE
1639

STERILE	EO
---------	----




4801141

Drape with hole
and adhesive
50x60cm




4803226

Folded drape
50x50cm




4803237


Slit drape with hole
60x60cm



4803239

Folded drape
60x60cm






STERILE EO




4802152

PP graduated cup,
round
120ml




4802153

PP graduated cup,
round
250ml




4802154

PP graduated cup,
round
500ml



4802155


PP graduated cup,
round
1000ml




4909976

Large measuring set

- 1x white tray
(258x153x45mm)
- 1x 500ml jug
- 1x 300ml cup
- 3x 60ml cups
- **Wrapped in 1
drape 60x60cm**






STERILE EO




4802122

Anatomical tweezers
13cm




4802123

Pean forceps 13 cm




8006001

Sponge forceps
24cm




8001734


Staple remover



8001740

Curette, 2 heads
14cm






STERILE EO




8001801

Scalpel N°15




8001806


Scalpel N°11



8001807

Scalpel N°10





SCISSORS


CE
1639

STERILE EO




8001749

Scissors 14cm
Blunt/blunt




8001712

Scissors 14cm
Sharp/blunt




8001716

Scissors 14cm
Sharp/sharp





8001817

Lister bandage scissors
14cm




8001820


Lister bandage scissors
18cm



8001720

Iris scissors
11.5cm





SCISSORS


CE
1639

STERILE EO




8001753

Curved Iris scissors
11.5cm




8001747

Fine Joseph scissors
15cm




8001748

Curved Jameson scissors
15cm




8001744

Curved Stevens scissors
11.5cm



8001757

Curved Metzenbaum scissors
18cm






HOOK - CURETTES


CE
1639

STERILE **EO**





4901087

Gillies hook
18cm




8001852

Curette 14cm
Head 3mm





8001750

Curette 14cm
Head 4mm





8001850

Curette 14cm
Head 7mm




8001765


ENT Curette



8001707

Blackhead remover curette
14cm





RETRACTORS – SPECULUM


CE
1639

STERILE EO





8001758

Holzheimer retractor
Blunt, 10cm




8001791 | 8001795

Senn-Miller retractor
Fine, 3 teeth




8001792


Weitlaner retractor
Pointed 10.5cm



8001623

Killian Speculum
35mm





FORCEPS


CE
1639

STERILE EO




8001705

Mosquito forceps
12.5cm




8001702

Curved Mosquito forceps
12.5cm




8001703

Hemostatic forceps
14-14.5cm




8001041

Hemostatic forceps
16cm




8001732


Kocher forceps
14cm



8001030

Kocher forceps
16cm





FORCEPS


CE
1639

STERILE EO




8001759

Crocodile forceps
14cm




8001600

Crocodile forceps Hartmann
4mm




8001601

Crocodile forceps Hartmann
6mm




8001602

Crocodile forceps Hartmann
8mm




8001616


Magill forceps
20cm



8001617

Magill forceps
25cm





FORCEPS - TWEEZERS


CE
1639

STERILE **EO**




8001620

Bayonet forceps
14cm





8001621

Bayonet forceps
16.5cm




8001619

Bayonet forceps
12.5cm




8001724


Adson tweezers 12cm
with claws




8001745

Adson tweezers 12cm
with fine claws







FORCEPS - TWEEZERS




8001742

Adson tweezers 12cm
Blunt





8001755

Adson tweezers 12cm
Blunt, fine





8001754

Feilchenfeld splinter tweezers
9cm




8001723

Surgical tweezers
14cm




8001721


Anatomical tweezers 13cm
Fine



8001731

Anatomical tweezers
14-14.5cm





**NEEDLE HOLDERS –
BLADES**


CE
1639

STERILE EO




8001718

Needle holder
14cm





8001719

Needle holder 14cm
Fine




8001502

Converse needle holder
(Tungsten clamp)
13cm




8001706


Beaver blade holder



8001710

Blade holder
Scalpel N°3






GYNAECOLOGY




4801060

Cusco vaginal speculum
XS – 75x17mm




4801061

Cusco vaginal speculum
S – 80x32mm




4801062

Cusco vaginal speculum
M – 100x25mm



4801063

Cusco vaginal speculum
L – 105x27mm






GYNAECOLOGY



CE
1639

STERILE EO





8001797

Anatomical forceps
24cm




8001630

Foerster sponge forceps
25cm





8001777

Pelkmann forceps
25cm




8001778


Curved Metzenbaum
Scissors, fine, 25cm



8001779


Pozzi forceps
25cm






WASHING PROCESS OF METAL INSTRUMENTS

Client satisfaction and product quality remain our top priorities. This is the reason why we have implemented a washing process and invested in washer-disinfectors to ensure that our metal instruments meet our machining residues and microbiological contamination acceptance criteria.





High-frequency surgery

Bipolar forceps are coagulation surgical instruments which :


- ✓ Minimise surgery trauma,
- ✓ Allow innovative solutions (no foreign body in the wound),
- ✓ Provide minimal and short interventions in the interests of patient's comfort,
- ✓ Reduce surgery costs.



8001737

Adson forceps
Double articulated connection


Tip width 0,5mm
Forceps length 12cm
Cable length 3m



8001818

Adson Forceps
Erbe connection


Tip width 0,5mm
Forceps length 12cm
Cable length 3m



8001819

Jeweler fine forceps
Double articulated connection

Tip width 0,5mm
Forceps length 11.5cm
Cable length 3m





High-frequency surgery

Bipolar forceps are coagulation surgical instruments which :


- ✓ Minimise surgery trauma,
- ✓ Allow innovative solutions (no foreign body in the wound),
- ✓ Provide minimal and short interventions in the interests of patient's comfort,
- ✓ Reduce surgery costs.



8001767

Straight-angled non-stick forceps
Double articulated connection


Tip width 0,5mm
Forceps length 12cm
Cable length 3m



8001768

Straight-angled non-stick forceps
Double articulated connection


Tip width 1mm
Forceps length 15.5cm
Cable length 3m




8001793

Straight-angled non-stick Jeweler
fine forceps
Double articulated connection

Tip width 0,5mm
Forceps length 10.8cm
Cable length 3m





High-frequency surgery

Bipolar forceps are coagulation surgical instruments which :


- ✓ Minimise surgery trauma,
- ✓ Allow innovative solutions (no foreign body in the wound),
- ✓ Provide minimal and short interventions in the interests of patient's comfort,
- ✓ Reduce surgery costs.



8001936

Bayonet shape non-stick forceps
Double articulated connection


Tip width 1mm
Forceps length 15.5cm
Cable length 3m





8001937

Bayonet shape non-stick forceps
Double articulated connection

Tip width 1mm
Forceps length 18cm
Cable length 3m





Stainless steel


CE
1639

STERILE EO




8001726

Olive tip cannula
Laser-welded olive
60mm / 19G




8001727

Olive tip cannula
Laser-welded olive
100mm / 18G




8001738

Olive tip cannula
Laser-welded olive
80mm / 17G




8001752

Olive tip cannula
Laser-welded olive
80mm / 14.5G




8001761

Olive tip cannula
Laser-welded olive
80mm / 18G




8001762

Olive tip cannula
Laser-welded olive
40mm / 17G




8001763

Olive tip cannula
Laser-welded olive
35mm / 21G




8001788

Olive tip cannula
Laser-welded olive
80mm / 21G




8001736

Olive tip cannula
Laser-welded olive
100mm / 18G
Double packaging




8001760

Olive tip cannula
Laser-welded olive
100mm / 16G




8001789



Olive tip cannula
Laser-welded olive
100mm / 19G



8001798

Olive tip cannula
Laser-welded olive
100mm / 13G





Stainless steel


CE
1639

STERILE EO




8001841

Olive tip cannula
Laser-welded olive
120mm / 13G




8001842

Olive tip cannula
Laser-welded olive
120mm / 14G




8001843

Olive tip cannula
Laser-welded olive
120mm / 15G




8001844

Olive tip cannula
Laser-welded olive
120mm / 17G




8001845



Curved olive tip cannula
Laser-welded olive
55mm / 17G



8001988

Olive tip cannula
Laser-welded olive
150mm / 18G





Stainless steel


CE
1639

STERILE EO




8001741

Blunt tip cannula
80mm / 16G




8001770

Curved blunt tip cannula
80mm / 19G




8001771

Curved blunt tip cannula
100mm / 23G




8001772

Curved blunt tip cannula
90mm / 17G




8001773

Curved blunt tip cannula
90mm / 19G




8001774

Curved blunt tip cannula
100mm / 21G




8001780

Curved blunt tip cannula
80mm / 18G




8001781

Curved blunt tip cannula
80mm / 16G




8001783

Curved blunt tip cannula
80mm / 14G




8001785

Curved blunt tip cannula
80mm / 13G




8001786

Blunt tip cannula
80mm / 13G



8001787

Blunt tip cannula
100mm / 11G






Plastic




8001730


Olive tip cannula
70mm



8001756

Olive tip cannula
100mm





COLLECTION AND RECYCLING OF SINGLE-USE METAL INSTRUMENTS


Applimed SA can supply its customers with more than a hundred single-use metal instruments. Since we are concerned about the impact on the environment, we provide our customers with a free collection and recycling service for our discarded metal instruments (after use).

The process of collection and recycling is described below:

-  Free supply of single-use 30-litre, 12-litre and 7-litre containers.
(these are specially designed for medical care waste resulting from a medical activity according to UN/OMoD 180101/180102 standards which comply with federal oMoD directives).
-  After use, the metal instruments for disposal are placed in the single-use container through a resealable cap. It is therefore no longer necessary to decontaminate them before placing them in the container.
-  Once full *(do not exceed the load indicated on the container)*, the container is hermetically sealed by the customer *(according to the instructions of use)*.
-  The company Rudosa collects the full containers which are hermetically sealed and exchanges them for new ones.
-  Full containers are collected and then transported according to ADR/SDR standards in an approved facility for incineration and/or ecological recovery.
-  This installation allows the containers to be introduced directly into the incinerator *(without opening them)* thus avoiding any risk of contamination.
-  All plastic components and residues on the instruments are burned and the metal completely decontaminated by the incineration is recovered by a physical system *(courant de Foucault)*.
-  The metal instruments are then sent to a different Swiss company specialised in metal recycling. They are combined with other metals in order to be recycled for technical uses.



This service is completely free and no commercial transaction applies. The sole purpose of this service is to facilitate the disposal of metals.




APPLIMED SA

SURGICAL FACE MASKS



9501001

Surgical face mask -
type IIR
CHILDREN



9501005

Surgical face mask -
type IIR
ADULT






GOWNS – OVERALLS




9501003

White gown
Size XL




9501007

Green protective gown
Size XL




9501008

Green plastified gown
Knitted cuffs
Size XL




9501053

Blue patient gown
Velcro opening at the front
Size M




9501004

Short white overall
140x77cm




9501022

Opaque blue overall
113x77cm



9501054

Blue patient gown
Velcro opening at the front
Size L






ANTI-SLIP SHOE COVER CAP PROTECTIVE CLOTHING SET



9501009


Blue anti-slip shoe cover
18x40cm



9508006


Protective clothing Set


- 1 dark blue overall
- 1 dark blue underwear
 - 1 cap
- 1 pair of blue anti-slip shoe covers



9501015

Cap
One size fits all






9507000

Bed protective set


- 1 mattress cover
220x96cm
- 1 bed sheet
220x160cm
- 1 pillowcase
50x70cm



9507021

Mattress cover


with elastic
220x96cm



9507022

Non woven


bed sheet
220x160cm




9507023

Pillowcase

60x60cm







Set composition

- 1 Blue plastic tweezers 13cm
 - 4 Gauze balls Ø 2.5-3cm
- Final packaging: rigid blister pack with 3 compartments/paper

Reference
4802116

STERILE	EO
---------	----






Set composition


- 3 Gauze balls Ø 2.5-3cm
 - 2 8-ply woven compresses, 10x10cm
 - 1 Iris scissors 11.5cm
 - 1 Anatomical tweezers 14-14.5cm
- Final packaging: rigid blister pack with 3 compartments/paper

Reference
4908808

STERILE	EO
---------	----




CE
1639




Set composition

- 6 Gauze balls Ø 2.5-3cm
 - 2 Non-woven compresses 7.5x7.5cm
 - 1 Drape 50x50cm
 - 1 Iris scissors 11.5cm
 - 1 Adson tweezers, blunt 12cm
- Final packaging: rigid blister pack with 3 compartments/paper

Reference
4909683






Set composition


- 1 Cup 250ml
 - 2 Anatomical tweezers 14-14.5cm
- Wrapped in a 60x60cm drape**
Final packaging: soft film peel pack/paper

Reference
4909798

STERILE	EO
---------	----



CE
1639




Set composition



- 1 Tray with 2 compartments, 180x95x25mm
- 1 Green plastic tweezers 13cm
- 1 Tray 60ml
- 5 Gauze balls Ø 2.5-3cm

Wrapped in a 50x50cm drape

Final packaging: soft film peelpack/paper

**Reference
4909927**






Set composition

- 2 Non-woven compresses 5x5cm
 - 4 Wooden applicators, big head, American shape
 - 1 Drape 50x50cm
- Final packaging: rigid blister pack with 3 compartments/paper

Reference
4909950






Set composition

- 1 Blue tray 210x140x45mm
 - 10 Gauze balls Ø 5-7cm
 - 1 Foerster sponge forceps 25cm
- Wrapped in a 80x80cm drape**
Final packaging: soft film peelpack/paper

Reference
4901045

STERILE	EO
---------	----



CE
1639




Set composition


- 1 White tray with 2 compartments, 750ml
- 6 Gauze balls Ø 2.5-3cm
- 1 Blue plastic tweezers 13cm
- 1 Drape 75x90cm

Wrapped in a 80x80cm drape

Final packaging: soft film peelpack/paper

**Reference
4901090**






Set composition


- 1 Cup 300ml
 - 5 Gauze balls Ø 5-7cm
 - 1 Foerster sponge forceps 25cm
 - 1 Drape 50x50cm
- Final packaging: soft film peelpack/paper

Reference
4909675

STERILE	EO
---------	----



CE
1639




Set composition

- 1 Sponge forceps 24 cm
 - 1 Cup 500ml
 - 6 Gauze balls **with X-ray thread**, Ø 5-7cm
- Final packaging: soft film peelpack/paper

Reference
4909677

STERILE	EO
---------	----



CE
1639




Set composition

- 1 Tray 125ml
 - 4 Gauze balls Ø 4-4.5cm
 - 1 Sponge forceps 18-19cm
- Final packaging: soft film peelpack/paper

Reference
4908805

STERILE	EO
---------	----




CE
1639



APPLIMED SA

SKIN DISINFECTION/ PREPARATION SET



Set composition


- 10 Non-woven balls Ø 4.5-5cm
- 1 Plastic sponge forceps 24cm

2 traceability labels


Primary packaging: rigid blister pack with 3 compartments/paper

Final packaging: soft film peelpack/paper or Steriking pouch

Reference
4901895




CE
1639



Set composition

- 3 Gauze balls Ø 2.5-3cm
 - 3 Woven compresses 7.5x7.5cm
 - 1 Blue plastic tweezers 13cm
- Final packaging: rigid blister pack with 3 compartments/paper

Reference
4801115






Set composition

- 1 Tray
 - 8 Gauze balls Ø 4cm
 - 2 Non-woven compresses 5x5cm
 - 1 Blue plastic tweezers 13cm
 - 1 Green plastic tweezers 13cm
- Wrapped in a 50x50cm drape**
Final packaging: soft film peelpack/paper

Reference
4901028

STERILE	EO
---------	----






Set composition


- 1 Tray
 - 5 Gauze balls **with X-ray thread** Ø 3cm
 - 3 **12-ply** woven compresses, 7.5x10cm
 - 3 **8-ply** woven compresses, 5x5cm
 - 3 **8-ply** woven compresses **with incision**, 5x5cm
 - 1 Blue plastic tweezers 13cm
 - 1 Cup 60ml
 - 1 Drape 50x50cm
- Wrapped in a 50x50cm drape**
Final packaging: soft film peelpack/paper

Reference
4809901

STERILE	EO
---------	----



CE
1639




Set composition


- 5 Gauze balls Ø 2.5-3cm
 - 5 Non-woven compresses 7.5x7.5cm
 - 1 Iris scissors 11.5cm
 - 1 Adson tweezers, blunt 12cm
 - 1 Blue plastic tweezers 13cm
 - 1 Green plastic tweezers 13cm
 - 1 Drape 50x50cm
- Final packaging: rigid blister pack with 3 compartments/paper

Reference
4901096

STERILE	EO
---------	----



CE
1639




Set composition



- 8 Non-woven balls Ø 3-3.5cm
 - 12 4-ply non-woven compresses 7.5x7.5cm
 - 1 Scissors sharp/blunt 14cm
 - 1 Anatomical tweezers 14-14.5cm
 - 1 Adson tweezers, blunt 12cm
 - 1 Drape 50x50cm
- Final packaging: rigid blister pack with 3 compartments/paper

Reference
4909083

STERILE	EO
---------	----




CE
1639




Set composition

- 2 8-ply woven compresses, 5x5cm
 - 5 8-ply woven compresses, 10x10cm
- Final packaging: soft film peelpack/paper

Reference
4909745



CE
1639




Set composition


- 3 Gauze balls Ø 2.5-3cm
 - 3 Non-woven compresses, 7.5x7.5cm
 - 1 Adson tweezers fine blunt 12cm
 - 1 Iris Scissors 11.5 cm
- Final packaging: : rigid blister pack with 3 compartments/paper

Reference
4909808

STERILE	EO
---------	----



CE
1639




Set composition


- 4 Gauze balls Ø 2.5-3cm
 - 4 4-ply non-woven compresses, 10x20cm
 - 1 Scissors, sharp/blunt 14cm
 - 1 Anatomical tweezers 14-14.5cm
- Final packaging: soft film peelpack/paper

Reference
4909857

STERILE	EO
---------	----



CE
1639




Set composition


- 1 Rigid blister pack with 2 compartments, 130x200mm
 - 10 Gauze balls Ø 2.5-3cm
 - 5 Non-woven compresses 7.5x7.5cm
 - 1 Scissors, sharp/blunt 14cm
 - 1 Iris scissors 11.5cm
 - 1 Anatomical tweezers 14cm
 - 1 Adson tweezers blunt 12cm
 - 1 Scalpel n°15
 - 1 Folded drape 60x60cm
- Wrapped in a 50x50cm drape**
Final packaging: soft film peelpack/paper

Reference
4909841

STERILE	EO
---------	----



CE
1639





Set composition

- 1 Blue plastic tweezers 13cm
 - 1 Green plastic tweezers 13cm
 - 5 Gauze balls Ø 2.5-3cm
 - 5 Non-woven compresses 7.5x7.5cm
 - 1 Folded drape 50x40cm
- Final packaging: rigid blister pack with 3 compartments/paper

Reference
4909920

STERILE	EO
---------	----








Set composition

- 1 Green plastic tweezers 13cm
 - 1 Blue plastic tweezers 13cm
 - 6 Gauze balls Ø 2.5-3cm
 - 2 Non-woven compresses 7.5x7.5cm
 - 1 Folded drape 50x40cm
- Final packaging: rigid blister pack with 3 compartments/paper

Reference
4909942

STERILE	EO
---------	----






Set composition

- 1 Tray with 2 compartments 180x95mm
 - 3 Non-woven balls Ø 3-4cm
 - 1 Blue plastic tweezers 13cm
 - 1 Green plastic tweezers 13cm
 - 4 Non-woven compresses 7.5x7.5cm
 - 1 Drape 50x50cm
- Wrapped in a 50x50cm drape**
Final packaging: soft film peelpack/paper

Reference
4909945

STERILE	EO
---------	----






Set composition


- 1 Drape with hole and adhesive 50x60cm
 - 1 Drape 75x80cm
 - 10 8-ply woven compresses, 10x10cm
 - 10 Gauze balls Ø 2.5-3cm
 - 5 8-ply woven compresses, 5x5cm
 - 1 2-piece Luer syringe, 2ml
 - 1 2-piece Luer syringe, 10ml
 - 1 Pink needle 38mm 18G
 - 1 Black needle 30mm 22G
 - 2 Tray 125ml
- Wrapped in a 80x80cm drape**
Final packaging: soft film peelpack/paper

Reference
4909979

STERILE	EO
---------	----




CE
1639




Set composition

- 3 Woven balls **(sterile)** Ø 2.5-3cm (4102008)
 - 1 Pink needle **(sterile)** 40mm, 18G
 - 1 Vasofix IV cannula pink **(sterile)** 20G
 - 1 Tegaderm **(sterile)** 7x8.5cm
 - 1 Seringue, Luer 3 pièces **(sterile)** 10ml
 - 1 3-way stopcock with 10cm extension **(sterile)**
 - 1 Hemostatic dressing Tip-Stop **(sterile)**
 - 2 Nitrile gloves with hem, size M
 - 2 Drape 50x50cm
 - 1 Elastic band 3x1.2mm
- Final packaging: soft film peelpack/paper

Reference
9006195






Set composition


- 1 Kidney-shaped APET bowl 252x150x40mm
- 1 Vasofix IV cannula pink **(sterile)** 20G
- 1 3-way stopcock with 10cm extension **(sterile)**
- 5 Woven balls **(sterile)** Ø 2.5-3cm (4102002)
- 1 Applica IV **(sterile)** 80x60 mm
- 1 0.9% NaCl syringe **(sterile)** 10ml
- 2 Nitrile gloves with hem, size M
- 1 Drape 50x50cm

Wrapped in a 60x60cm drape

Final packaging: soft film peelpack/paper

Reference
9009186






Set composition



- 1 Kidney-shaped APET bowl 252x150x40mm
- 1 Vasofix IV cannula green **(sterile)** 18G
- 1 3-way stopcock with 10cm extension **(sterile)**
- 5 Woven balls **(sterile)** Ø 2.5-3cm (4102002)
- 1 Applica IV **(sterile)** 80x60 mm
- 1 0.9% NaCl syringe **(sterile)** 10ml
- 2 Nitrile gloves with hem, size M
- 1 Drape 50x50cm

Wrapped in a 60x60cm drape

Final packaging: soft film peelpack/paper

Reference
9009187






Set composition

- 1 Blue tray 210x140x45mm
 - 1 Needle holder 14cm
 - 1 Kocher forceps 14cm
 - 1 Baby Metzenbaum scissors 11.5cm
 - 1 Adson tweezers with fine claws 12cm
 - 1 Anatomical tweezers 16cm
 - 2 Cup 120ml
 - 5 Non-woven balls $\varnothing$ 4.5cm
 - 1 Hemostatic forceps 14-14.5cm
 - 1 Olive tip cannula 100mm, 18G
 - 1 Drape with hole and adhesive 75x90cm
- Wrapped in a 80x80cm drape**
Final packaging: soft film peel pack/paper

Reference 4901047



1639






Set composition

- 1 White tray with 2 compartments 750ml
 - 6 Gauze balls Ø 2.5-3cm
 - 5 Non-woven compresses 5x5cm
 - 3 Non-woven compresses 10x10cm
 - 1 Iris scissors 11.5cm
 - 1 Adson tweezers blunt 12cm
 - 1 Anatomical tweezers 14-14.5cm
 - 1 Needle holder 14cm
 - 1 Pink needle 38mm 18G
 - 1 Orange needle 25mm 25G
 - 1 2-piece Luer syringe 5ml
 - 2 Nitrile gloves with hem, size M
 - 1 Drape with hole and adhesive 50x60cm
- Wrapped in a 60x60cm drape**
Final packaging: soft film peel pack/paper

Reference
4901089

STERILE EO






Set composition

- 5 4-ply non-woven compresses, 7.5x7.5cm
 - 1 Iris scissors 11.5cm
 - 1 Adson tweezers with claws 12cm
 - 1 Needle holder 14cm
 - 1 Folded drape 50x50cm
 - 1 Folded drape with hole and adhesive 50x60cm
- Final packaging: rigid blister pack with 3 compartments/paper

Reference
4903236

STERILE	EO
---------	----



CE
1639




Set composition


- 1 White tray with 2 compartments 750ml
 - 8 Gauze balls Ø 2.5-3cm
 - 5 Non-woven compresses 5x5cm
 - 5 Non-woven compresses 10x10cm
 - 1 Curved Stevens scissors 11.5cm
 - 1 Adson tweezers with fine claws 12cm
 - 1 Fine anatomical tweezers, 13cm
 - 1 Converse needle holder 12.5cm
 - 1 Scalpel n°15
 - 1 Drape with hole and adhesive 50x60cm
- Wrapped in a 60x60cm drape**
Final packaging: soft film peelpack/paper

Reference
4909850

STERILE	EO
---------	----




CE
1639





Set composition

- 1 Rigid blister pack with 2 compartments, 130x200mm
 - 6 Gauze balls Ø 2.5-3cm
 - 5 Non-woven compresses 5x5cm
 - 1 Metzenbaum scissors curved 15cm
 - 1 Adson tweezers with claws 12cm
 - 1 Adson tweezers blunt 12cm
 - 1 Anatomical tweezers 14cm
 - 1 Needle holder 14cm
 - 1 Drape with hole 75x80cm
 - 1 **Wrapped in a 60x60cm drape**
- Final packaging: soft film peel pack/paper

Reference
4909889






Set composition

- 1 Rigid blister pack with 2 compartments, 130x200mm
 - 6 Gauze balls Ø 2.5-3cm
 - 5 Non-woven compresses 5x5cm
 - 1 Needle holder 14cm
 - 1 Iris scissors 11.5cm
 - 1 Adson tweezers blunt 12cm
 - 1 Adson tweezers with claws 12cm
 - 2 Nitrile gloves with hem, size M
 - 1 Drape with hole and adhesive 50x60cm
 - 1 Orange needle 25mm, 25G
 - 1 Pink needle 38mm, 18G
 - 1 3-piece Luer syringe, 5ml
 - 1 Mask
 - 1 Cellulose mattress protector 40x60cm
- Wrapped in a 60x60cm drape**
Final packaging: soft film peelpack/paper

Reference
4909961






Set composition


- 1 Scissors 14cm sharp/sharp
 - 1 Anatomical tweezers 14-14.5cm
 - 1 Adson tweezers with claws 12cm
 - 1 Needle holder 14cm
- Final packaging: rigid blister pack with 3 compartments/paper

Reference
4909962

STERILE	EO
---------	----



CE
1639



Set composition


- 1 Rigid blister pack with 2 compartments, 130x200mm
- 6 Gauze balls Ø 2.5-3cm
- 5 Non-woven compresses 5x5cm
- 1 Needle holder 14cm
- 1 Iris scissors 11.5cm
- 1 Adson tweezers with claws 12cm
- 1 Adson tweezers blunt 12cm
- 1 Drape with hole and adhesive 50x60cm

Wrapped in a 50x50cm drape



Final packaging: soft film peel pack/paper

Reference
4909968

STERILE	EO
---------	----



CE
1639




Set composition

- 1 Rigid blister pack with 2 compartments, 130x200mm
- 6 Gauze balls Ø 2.5-3cm
- 5 Non-woven compresses 5x5cm
- 1 Needle holder 14cm
- 1 Iris scissors 11.5cm
- 1 Adson tweezers with fine claws 12cm
- 1 Pean forceps 13cm
- 1 Adhesive drape with hole 50x60cm



Wrapped in a 80x80cm drape

Final packaging: soft film peelpack/paper

Reference
4909969



CE
1639




Set composition



- 1 Rigid blister pack with 2 compartments, 130x200mm
 - 8 Gauze balls Ø 2.5-3cm
 - 10 Non-woven compresses 5x5cm
 - 1 Derf needle holder 12.5cm
 - 1 Iris scissors 11.5cm
 - 1 Fine anatomical tweezers 13cm
 - 1 Adson tweezers with fine claws 12cm
 - 1 Curved Stevens scissors 11.5cm
 - 1 Pink needle 38mm, 18G
 - 1 Orange needle 25mm, 25G
 - 1 3-piece Luer syringe, 5ml
 - 1 Scalpel n°15
 - 1 Drape with hole and adhesive 50x60cm
- Wrapped in a 60x60cm drape**
Final packaging: soft film peel pack/paper

Reference
4909971

STERILE	EO
---------	----



CE
1639




Set composition

- 1 White tray with 2 compartments 750ml
 - 4 Gauze balls Ø 2.5-3cm
 - 4 Non-woven compresses 5x5cm
 - 4 Non-woven compresses 10x10cm
 - 1 Iris scissors 11.5cm
 - 1 Adson tweezers with fine claws 12cm
 - 1 Mosquito forceps, curved 12.5cm
 - 1 Converse needle holder 12.5cm
 - 1 Drape with hole and adhesive 50x60cm
- Wrapped in a 60x60cm drape**
Final packaging: soft film peelpack/paper

Reference
4909993

STERILE	EO
---------	----




CE
1639





Set composition

- 1 Blue tray 210x140x45mm
- 1 Tray 125ml
- 1 Cup 60ml
- 1 Red needle 40mm, 18G
- 1 Orange needle 25mm, 25G
- 1 2-piece Luer syringe, 5ml
- 1 3-piece Luer syringe, 10ml
- 5 Gauze balls Ø 2.5-3cm
- 5 8-ply woven compresses, 5x5cm
- 4 12-ply woven compresses, 7.5x7.5cm
- 1 Drape with hole and adhesive 75x90cm
- 1 Cellulose mattress protector 40x60cm
- 1 Converse needle holder 12.5cm
- 1 Mosquito forceps, curved 12.5cm
- 1 Hemostatic forceps 14-14.5cm
- 1 Adson tweezers with fine claws 12cm
- 1 Jameson scissors, fine 15cm
- 1 Temporary sharps holder
- Wrapped in a 75x75cm Dextex drape**
- Sterile packaging – 1st layer : soft film peelpack/paper**
- 2 Latex gloves «Vasco OP Sensitive», size 7.5 **sterile**
- Final packaging : Steriking bag**

Reference
9009203



CE
1639




Set composition


- 6 Gauze balls Ø 2.5-3cm
 - 2 4-ply non-woven compresses, 7.5x7.5cm
 - 1 Iris scissors 11.5cm
 - 1 Adson tweezers blunt 12cm
 - 1 Drape 50x50cm
- Final packaging: rigid blister pack with 3 compartments/paper

Reference
4909883

STERILE	EO
---------	----



CE
1639




Set composition


- 4 Gauze balls Ø 2.5-3cm
 - 1 Adson tweezers blunt 12cm
 - 1 Scalpel n°11
- Final packaging: rigid blister pack with 3 compartments/paper

Reference
4909936

STERILE	EO
---------	----



CE
1639



APPLIED SA


SUTURE REMOVAL SET



Set composition

- 1 Iris scissors 11.5cm
 - 1 Adson tweezers blunt 12cm
- Final packaging: soft film peelpack/paper

**Reference
4909996**



CE
1639



Set composition


- 1 Kidney-shaped APET bowl 252x150x40mm
- 1 Tray 125ml
- 6 Gauze balls Ø 2.5-3cm
- 1 Blue tweezers 13cm
- 1 Drape with hole and slit, 60x60cm
- 3 Nitrile gloves with hem, size L

Wrapped in a 60x60cm drape

Final packaging: soft film peelpack/paper

**Reference
4908814**

STERILE	EO
---------	----



CE
1639



Set composition


- 1 Kidney-shaped APET bowl 252x150x40mm
- 6 Gauze balls Ø 2.5-3cm
- 1 Tray 125ml
- 1 Blue plastic tweezers 13cm
- 1 Drape with hole and slit, 60x60cm
- 1 Drape 60x60cm

Wrapped in a 60x60cm drape

Final packaging: soft film peelpack/paper

**Reference
4908823**

STERILE	EO
---------	----



CE
1639



Set composition


- 1 Kidney-shaped APET bowl 252x150x40mm
- 6 Gauze balls Ø 2.5-3cm
- 1 Blue plastic tweezers 13cm
- 1 Tray 125ml
- 2 Nitrile gloves with hem, size M
- 1 Drape with hole and slit 60x60cm

Wrapped in a 60x60cm drape


Final packaging: soft film peelpack/paper

**Reference
4809913**

STERILE	EO
---------	----



CE
1639




Set composition



- 1 Tray 125ml
 - 8 8-ply woven compresses, 5x5cm
- Final packaging: soft film peelpack/paper

**Reference
4909906**

STERILE	EO
---------	----



CE
1639




Set composition


- 1 Rigid blister pack with 2 compartments, 130x200mm
 - 5 Woven balls Ø 2.5-3cm
 - 5 4-ply non-woven compresses, 5x5cm
 - 1 Scalpel n°3 blade holder
 - 1 curved Iris scissors 11.5cm
 - 1 Converse needle holder 12.5cm
 - 1 Adson tweezers fine blunt 12cm
 - 1 Adson tweezers with fine claws 12cm
 - 1 3-piece LL syringe, 3ml
 - 1 Pink needle 38mm, 18G
 - 1 Yellow needle 13mm, 30G
 - 1 Drape with hole and adhesive 50x60cm
- Wrapped in a 50x50cm drape**
Final packaging: soft film peelpack/paper

Reference
4909991

STERILE	EO
---------	----



CE
1639




Set composition



- 1 Blue tray 210x140x45mm
 - 2 Cup 60ml
 - 10 Non-woven compresses 7.5x7.5cm
 - 5 Gauze balls Ø 2.5-3cm
 - 1 Mosquito forceps, curved 12.5cm
 - 1 Adson tweezers fine blunt 12cm
 - 1 Adson tweezers with fine claws 12cm
 - 1 Iris scissors 11.5cm
 - 1 Converse needle holder 12.5cm
 - 1 Baby Metzenbaum scissors 11.5cm
 - 1 Scalpel n°15
 - 1 Blue needle 40mm, 23G
 - 1 Orange needle 25mm, 25G
 - 1 3-piece LL syringe, 5ml
 - 1 Drape with hole and adhesive 75x90cm
- Wrapped in a 80x80cm drape**
Final packaging: soft film peel pack/paper

Reference
4908826

STERILE	EO
---------	----



CE
1639





Set composition

- 1 White tray with 2 compartments, 750ml
 - 6 Gauze balls Ø 2.5-3cm
 - 5 Non-woven compresses 5x5cm
 - 5 Non-woven compresses 10x10cm
 - 1 Converse needle holder 12.5cm
 - 1 Metzenbaum scissors curved 15cm
 - 1 Iris scissors 11.5cm
 - 1 Adson tweezers blunt 12cm
 - 1 Adson tweezers with fine claws 12cm
 - 1 2-piece Luer syringe, 5ml
 - 1 Pink needle 38mm, 18G
 - 1 Blue needle 25mm, 23G
 - 1 Scalpel n°3 blade holder
 - 1 Drape with hole and adhesive 50x60cm
 - 1 Cellulose mattress protector 40x60cm
- Wrapped in a 80x80cm drape**
Final packaging: soft film peel pack/paper

Reference
4908828

STERILE	EO
---------	----






Set composition

- 1 Rigid blister pack with 2 compartments, 130x200mm
- 5 Gauze balls Ø 2.5-3cm
- 5 Non-woven compresses 5x5cm
- 1 Scalpel n°3 blade holder
- 1 Curved Stevens scissors 11.5cm
- 1 Converse needle holder 12.5cm
- 1 Adson tweezers with fine claws 12cm
- 1 Adson tweezers fine blunt 12cm
- 1 3-piece LL syringe, 3ml
- 1 Pink needle 38mm, 18G
- 1 Yellow needle 13mm, 30G
- 1 Drape with hole and adhesive 50x60cm

Wrapped in a 50x50cm drape

Final packaging: soft film peelpack/paper

**Reference
4909643**






Set composition


- 1 Rigid blister pack with 2 compartments, 130x200mm
- 8 Gauze balls Ø 2.5-3cm
- 1 Needle holder 14cm
- 1 Iris scissors 11.5cm
- 1 Scissors blunt 14cm
- 1 Curved Stevens scissors 11.5cm
- 1 Fine anatomical tweezers 13cm
- 1 Adson tweezers with fine claws 12cm
- 2 Gillies hook 18cm
- 2 Fine Senn-Miller retractors
- 1 Scalpel n°15
- 1 Scalpel n°10
- 1 3-piece Luer syringe, 5ml

Wrapped in a 60x60cm drape

Final packaging: soft film peelpack/paper

**Reference
4909093**







Set composition

- 1 White tray with 2 compartments, 750ml
 - 8 Gauze balls Ø 2.5-3cm
 - 6 4-ply non-woven compresses, 5x5cm
 - 7 4-ply non-woven compresses, 10x10cm
 - 1 Pad holding forceps 18-19cm
 - 1 Converse needle holder 12.5cm
 - 1 Iris scissors 11.5cm
 - 1 Scissors sharp/blunt 14cm
 - 1 Curved Mayo scissors 15cm
 - 1 Curved Stevens scissors 11.5cm
 - 1 Fine anatomical tweezers 13cm
 - 1 Adson tweezers with fine claws 12cm
 - 2 Gillies hook 18cm
 - 2 Fine Senn-Miller retractors
 - 2 Farabeuf retractors 12cm
 - 2 Rangell-Davis retractors 14cm
 - 2 Scalpel n°3 blade holders
 - 1 3-piece Luer syringe, 5ml
 - 1 Skin marker pen and 1 ruler
- Wrapped in a 60x60cm drape**
Final packaging: soft film peel pack/paper

Reference 4909094

STERILE	EO
---------	----





Set composition


- 1 Kidney-shaped APET bowl 252x150x40mm
- 1 Mosquito forceps 12.5cm
- 2 Curved Mosquito forceps 12.5cm
- 1 Needle holder 14cm
- 1 Cup 60ml
- 5 Gauze balls Ø 2.5-3cm
- 1 Pink needle 38mm, 18G
- 1 Brown needle 26mm, 26G
- 1 Black needle 30mm, 22G
- 1 3-piece Luer syringe, 5ml

Wrapped in a 60x60cm drape

Final packaging: soft film peelpack/paper

**Reference
4909816**

STERILE	EO
---------	----



CE
1639



Set composition


- 1 Kidney-shaped APET bowl 252x150x40mm
- 1 Blue needle 40mm, 23G
- 1 Orange needle 25 mm, 25G
- 1 Adson tweezers fine blunt 12cm
- 1 Double scraper tool 18.5cm
- 1 Converse needle holder 12.5cm
- 1 Kilner scissors 15cm
- 1 Mosquito forceps 12.5cm
- 1 2-piece Luer syringe, 10ml
- 2 Cups 60ml
- 1 Cup 250ml
- 8 Steel applicators, cotton head 150x0.9mm, Ø 2.2-2.7mm
- 1 Elastic band 3x1.2mm
- 5 Non-woven compresses 5x5cm
- 2 Non-woven compresses 10x10cm
- 1 Drape 75x90cm

Wrapped in a 80x80cm drape


Final packaging: soft film peel pack/paper

Reference
4908827

STERILE EO



CE
1639



Set composition


- 1 Kidney-shaped APET bowl 252x150x40mm
- 5 4-ply non-woven compresses, 5x5cm
- 2 4-ply non-woven compresses, 10x10cm
- 1 Nail scissors 15cm
- 1 Baby Metzenbaum scissors 11.5cm
- 1 Adson tweezers with fine claws 12cm
- 1 Kocher forceps 14cm
- 1 Hemostatic forceps 14-14.5cm
- 1 Needle holder 14cm
- 1 Scalpel n°3 blade holder
- 1 Double scraper tool
- 1 Curette 14cm, head 4mm
- 1 Cup 250ml
- 1 Elastic band
- 1 Drape 75x90cm

Wrapped in a 60x60cm drape



Final packaging: soft film peel pack/paper

Reference
4901097

STERILE EO



CE
1639



Set composition


- 1 Kidney-shaped APET bowl 252x150x40mm
- 1 Kocher forceps 14cm
- 1 Hemostatic forceps 14-14.5cm
- 1 Very fine surgical tweezers 14cm
- 1 Anatomical tweezers 14-14.5cm
- 1 Scissors sharp/blunt 14cm
- 1 Needle holder 18cm
- 10 Gauze balls Ø 2.5-3cm
- 10 12-ply woven compresses, 10x10cm
- 1 Drape 75x90cm

Wrapped in a 80x80cm drape


Final packaging: soft film peel pack/paper

**Reference
4909689**

STERILE	EO
---------	----



CE
1639




Set composition


- 2 Hemostatic forceps 14-14.5cm
 - 1 Mosquito forceps 12.5cm
 - 1 Scissors sharp/blunt 14cm
 - 1 Needle holder 18cm
 - 1 Anatomical tweezers 14-14.5cm
 - 1 Surgical tweezers 14cm
- Wrapped in a 60x60cm drape**
Final packaging: soft film peelpack/paper

Reference
4909684

STERILE	EO
---------	----



CE
1639




Set composition


- 1 Blue tray 210x140x45mm
 - 10 Gauze balls Ø 2.5-3cm
 - 1 Cup 250ml
 - 1 Pelkmann forceps 25cm
 - 1 Pozzi forceps 25cm
 - 1 Fine curved Metzenbaum scissors 25cm
 - 1 Cusco metal speculum 105x27mm
 - 1 Hysterometer (indicative scale) no measurement function
 - 2 Nitrile gloves with hem, size M
- Wrapped in a 60x60cm drape**
Final packaging: soft film peel pack/paper or Steriking pouch

Reference
4902214

STERILE	EO
---------	----



CE
1639




Set composition

- 1 Kidney-shaped APET bowl 252x150x40mm
- 10 16-ply woven compresses 10x10cm
- 1 Cellulose mattress protector 60x90cm
- 1 3-piece LL syringe **(sterile)** 20ml
- 1 Pink needle 38 mm **(sterile)** 18G
- 2 Surgical gown, size L **(sterile)**
- 1 Mini-Spike **(sterile)**


Wrapped in a 80x80cm drape

Final packaging: soft film peel pack/paper

**Reference
9001048**



CE
1639




Set composition

- 1 Kidney-shaped APET bowl 252x150x40mm
- 1 Tube KY 5G
- 10 16-ply woven compresses **(sterile)** 10x10cm
- 1 Cellulose mattress protector 60x90cm
- 1 3-piece LL syringe **(sterile)** 20ml
- 1 Pink needle 38mm **(sterile)** 18G
- 1 Blue needle 25mm **(sterile)** 23G
- 2 Green gowns, size L
- 1 Mini-Spike **(sterile)**
- 1 Sheet 80x210cm

Wrapped in a 80x80cm drape

Final packaging: soft film peel pack/paper


Reference
9009192





APPLIED SA

ENDOSCOPE CLEANING CARE SET



Set composition

- 1 Polythene bag 350x270x0.04mm
- 1 Cup 250ml
- 1 Endoscope tray protector
- 1 Cellulose mattress protector 60x90cm


Contained in a bag:

- 1 - Short brush 15.5cm
- 1 - Long brush 230cm
- 1 - Endoscope sponge

Wrapped in a 60x60cm drape

Final packaging: soft film peelpack/paper


**Reference
9009191**





APPLIED SA


ENDOSCOPE CLEANING CARE SET



Set composition

- 1 Short brush 15.5cm
 - 1 Long brush 230cm
 - 1 Endoscope sponge
- Final packaging: Minigrip bag

**Reference
9009184**





APPLIMED SA


CT SCANNER INJECTION SET



Set composition

- 1 Vasofix IV cannula pink **(sterile)** 20G
 - 3 Woven balls **(sterile)** Ø 2.5-3cm (4102008)
 - 1 3-way stopcock with 10cm extension **(sterile)**
 - 1 Applica IV **(sterile)** 80x60 mm
 - 1 Seringue, Luer 3 pièces **(sterile)** 10ml
 - 2 Nitrile gloves with hem, size M
 - 1 Elastic band 3x1.2mm
- Final packaging: soft film peelpack/paper

reference
9006066







Set composition

- 1 Blue tray 210x140x45mm
 - 1 Drape 75x80cm
 - 1 Drape with hole and adhesive 50x60cm
 - 2 Cups 250ml
 - 3 Non-woven balls Ø 3-4cm
 - 10 Non-woven compresses 5x5cm
 - 10 Non-woven compresses 7.5x7.5cm
 - 1 Pad holding forceps 24cm
 - 1 Curved Mosquito forceps 12.5cm
 - 2 Adson tweezers with fine claws 12cm
 - 1 Curved Stevens scissors 11.5cm
 - 1 Converse needle holder 12.5cm
 - 1 3-piece Luer syringe, 10ml
 - 1 Pink needle 38mm, 18G
 - 1 Orange needle 25mm, 25G
 - 1 Black needle 30mm, 22G
- Wrapped in a Dextex drape 90x90cm**
Final packaging: soft film peel pack/paper

**Reference
4803223**

STERILE	EO
---------	----







Set composition

- 1 Blue tray 210x140x45mm
- Individually packed inside a peelpack:
- 1 - Mosquito forceps 12.5cm
- 1 - Stevens scissors 12.5cm
- 1 - Scalpel n°3 blade holder
- 1 - Adson tweezers with fine claws 12cm
- 1 - Converse needle holder 12.5cm
- 1 - Adson tweezers fine blunt 12cm
- 1 - Gillies hook 18cm
- 1 traceability label**
- Wrapped in a Dextex drape 75x75cm**
- Final packaging: soft film peelpack/paper

**Reference
4909981**







Set composition

- 1 Blue tray 210x140x45mm
 - 5 4-ply non-woven compresses, 10x10cm
 - 1 Anatomical tweezers 14-14.5cm
 - 1 Mosquito forceps 12.5cm
 - 1 Olive tip cannula 60mm, 20G
 - 1 2-piece Luer syringe, 5ml
 - 1 2-piece Luer syringe, 10ml
 - 1 Cup 60ml
 - 1 Cup 120ml
 - 1 Scalpel n°11
 - 1 Drape with hole and adhesive 50x60cm
- Wrapped in a 60x60cm drape**
Final packaging: soft film peelpack/paper

Reference
4909914

STERILE EO






Set composition

- 1 White tray with 2 compartments, 750ml
- 5 Gauze balls Ø 2.5-3cm
- 1 Blue tweezers 13cm
- 6 8-ply woven compresses, 10x10cm
- 1 Folded drape 60x60cm

Wrapped in a 60x60cm drape

Final packaging: soft film peelpack/paper

**Reference
4909090**



CE
1639




Set composition

- 1 Rigid blister pack with 2 compartments, 130x200mm
- 6 Gauze balls Ø 2.5-3cm
- 5 Non-woven compresses 5x5cm
- 1 Converse needle holder 12.5cm
- 1 Iris scissors 11.5cm
- 1 Adson tweezers fine blunt 12cm
- 1 Adson tweezers with fine claws 12cm
- 1 Drape with hole and adhesive 50x60cm


Wrapped in a 80x80cm drape

Final packaging: soft film peelpack/paper

Reference
4909855



CE
1639




Set composition



- 2 Cups 120ml
 - 5 Gauze balls Ø 2.5-3cm
 - 5 Non-woven compresses 5x5cm
 - 2 Non-woven compresses 10x10cm
 - 1 Iris scissors 11.5cm
 - 1 Adson tweezers fine blunt 12cm
 - 1 Anatomical tweezers 14-14.5cm
 - 1 Scalpel n°3 blade holder
- Wrapped in a 45x45cm drape**
Final packaging: soft film peelpack/paper

Reference
4909862

STERILE	EO
---------	----



CE
1639




Set composition



- 1 Blue tray 210x140x45mm
 - 5 Non-woven compresses 10x20cm
 - 1 Cup 250ml
 - 3 Non-woven balls Ø 3-4cm
 - 1 3-piece Luer syringe, 20ml
 - 1 Pink needle 38mm, 18G
 - 1 Yellow needle 70mm, 20G
 - 1 Spencer Wells forceps 20cm
 - 1 Converse needle holder 12.5cm
 - 1 Scalpel n°11
 - 1 Very fine surgical tweezers 14cm
 - 1 Curved Metzenbaum scissors 18cm
 - 1 Pad holding forceps 24cm
 - 10 Non-woven compresses 7.5x7.5cm
- Wrapped in a Dextex drape 90x90cm**
Final packaging: soft film peelpack/paper

Reference
4803222

STERILE	EO
---------	----




CE
1639



Set composition

- 6 Gauze balls Ø 2.5-3cm
 - 2 Non-woven compresses 5x5cm
 - 2 Non-woven compresses 10x10cm
 - 1 Anatomical tweezers 14-14.5cm
 - 1 Cup 60ml
 - 1 Drape with hole and adhesive 50x60cm
- Wrapped in a 60x60cm drape**
1 traceability label
Final packaging: soft film peelpack/paper

Reference
4909917






Set composition

- 1 Rigid blister pack with 2 compartments, 130x200mm
- 4 Gauze balls Ø 2.5-3cm
- 5 Non-woven compresses 5x5cm
- 5 Non-woven compresses 10x10cm
- 1 Adson tweezers blunt 12cm
- 1 Drape with hole and adhesive 50x60cm

Wrapped in a 50x50cm drape

Final packaging: soft film peel pack/paper

**Reference
4909946**



CE
1639




Set composition


- 1 White tray with 2 compartments 750ml
 - 6 Gauze balls Ø 2.5-3cm
 - 1 Anatomical tweezers 14cm
 - 2 2-piece Luer syringes, 2ml
 - 1 2-piece Luer syringe, 5ml
 - 1 2-piece Luer syringe, 20ml
 - 1 Pink needle 38mm, 18G
 - 1 Green needle 40mm, 21G
 - 1 Blue needle 25mm, 23G
 - 1 Drape with hole and adhesive 50x60cm
- Wrapped in a 60x60cm drape**
Final packaging: soft film peel pack/paper

Reference
4909949

STERILE	EO
---------	----



CE
1639




Set composition


- 1 Blue tray 210x140x45mm
 - 5 Gauze balls Ø 2.5-3cm
 - 5 Non-woven compresses 10x20cm
 - 1 Scissors sharp/blunt 14cm
 - 2 Anatomical tweezers 14-14.5cm
 - 1 Hemostatic forceps 14-14.5cm
 - 1 Cup 250ml
- Wrapped in a 60x60cm drape**
Final packaging: soft film peel pack/paper

Reference
4909858

STERILE	EO
---------	----



CE
1639




Set composition

- 1 Blue tray 210x140x45mm
 - 1 Cup 250ml
 - 1 Cup 60ml
 - 1 Pad holding forceps 24cm
 - 4 12-ply woven compresses 10x20cm
 - 10 12-ply woven compresses, 10x10cm
 - 1 Iris scissors curved 11.5cm
 - 1 Adson tweezers blunt 12cm
 - 1 Adson tweezers with fine claws 12cm
 - 2 Mosquito forceps curved 12.5cm
 - 1 Fine needle holder 14cm
 - 1 Hemostatic forceps 14-14.5cm
 - 4 Gauze balls Ø 4-4.5cm
 - 1 Drape with hole and adhesive 75x90cm
- Wrapped in a 80x80cm drape**
Final packaging: soft film peel pack/paper


Reference
4909804

STERILE	EO
---------	----



CE
1639

APPLIMED SA



Address

Z.I. Pra de Plan 1
CH - 1618 Châtel-St-Denis

Contact

Reception : +41 21 948 92 74
Sales representative : +41 79 473 10 68
E-Mail : info@applimed.ch

Office hours

Monday – Thursday : 8am – 5pm
Friday : Closed
Saturday : Closed

www.applimed.ch