



High-frequency surgery

Bipolar forceps are coagulation surgical instruments which :


- ✓ Minimise surgery trauma,
- ✓ Allow innovative solutions (no foreign body in the wound),
- ✓ Provide minimal and short interventions in the interests of patient's comfort,
- ✓ Reduce surgery costs.



8001737

Adson forceps
Double articulated connection


Tip width 0,5mm
Forceps length 12cm
Cable length 3m



8001818

Adson Forceps
Erbe connection


Tip width 0,5mm
Forceps length 12cm
Cable length 3m



8001819

Jeweler fine forceps
Double articulated connection

Tip width 0,5mm
Forceps length 11.5cm
Cable length 3m





High-frequency surgery

Bipolar forceps are coagulation surgical instruments which :


- ✓ Minimise surgery trauma,
- ✓ Allow innovative solutions (no foreign body in the wound),
- ✓ Provide minimal and short interventions in the interests of patient's comfort,
- ✓ Reduce surgery costs.



8001767

Straight-angled non-stick forceps
Double articulated connection



Tip width 0,5mm
Forceps length 12cm
Cable length 3m



8001768

Straight-angled non-stick forceps
Double articulated connection


Tip width 1mm
Forceps length 15.5cm
Cable length 3m



8001793

Straight-angled non-stick Jeweler
fine forceps
Double articulated connection

Tip width 0,5mm
Forceps length 10.8cm
Cable length 3m





High-frequency surgery

Bipolar forceps are coagulation surgical instruments which :


- ✓ Minimise surgery trauma,
- ✓ Allow innovative solutions (no foreign body in the wound),
- ✓ Provide minimal and short interventions in the interests of patient's comfort,
- ✓ Reduce surgery costs.



8001936

Bayonet shape non-stick forceps
Double articulated connection

Tip width 1mm
Forceps length 15.5cm
Cable length 3m



8001937

Bayonet shape non-stick forceps
Double articulated connection

Tip width 1mm
Forceps length 18cm
Cable length 3m

